

## Bench-Top Solid Waste Container



- Made of chemically resistant HDPE for increased durability.
- Ideal containment system for chemically contaminated solid wastes such as, wet silica gel, radioisotopes, gloves, kimwipes, filter paper, wiping rags, pipettes, pipette bulbs, and other disposables.
- Prevents cross contamination of wastes within the laboratory.
- Has a small footprint and can easily sit in a corner of a fume hood or on the bench.
- Latching lid with gasket significantly reduces emissions of solvent-soaked solid wastes.

## Keep Your Laboratory Safe and Compliant


### Solid Waste Container

Item# SWC-002

Meets 40 CFR 264.172

Compatibility of Waste Containers:

"The owner or operator must use a container made of, or lined with, materials which will not react with, and are otherwise compatible with the hazardous waste to be stored, so that the ability of the container to contain the waste is not impaired."



Polyethylene Liners  
100 or 250 count  
Yellow or Red

### Solid Waste Container Bags

SWC-B100Y:	100 Yellow Liners	/	SWC-B100R:	100 Red Liners
SWC-B250Y:	250 Yellow Liners	/	SWC-B250Y:	250 Red Liners