

PolyPro RS-PB-100 & 200

Water Baths



Revolutionary Science

Manufacturer of Precision Laboratory Equipment

The Poly Pro Baths have unique polypropylene chambers. Due to the distinctive nature of polypropylene, not only is the chamber's surface resistant to corrosion and chemical attack, but the glossy finish is non-stick for an easy clean up. The digital display and controls make these units great for all general purpose laboratory applications! Two year warranty included.

RS-PB-100 (Single Chamber) specifications

Capacity:	5.5L
Net Weight:	2.72kg (6lbs)
Power Requirements:	100-240 VAC, 50-60 Hz, 150 W, 1.3 A
Dimensions:	
Exterior:	40cm wide x 25cm deep x 17cm tall 15.7" wide x 9.8" deep x 6.7" tall
Interior:	32cm wide x 17cm deep x 13.3cm tall 12.6" wide x 6.7" deep x 5.2" tall



RS-PB-200 (Dual Chamber) specifications

Capacity:	5.5L (each chamber)
Net Weight:	7.25kg (16lbs)
Power Requirements:	100-240 VAC, 50-60 Hz, 300 W, 2.6 A
Dimensions:	
Exterior:	56cm wide x 45.5cm deep x 14.7cm tall 22" wide x 17.9" deep x 5.8" tall
Interior:	32cm wide x 17cm deep x 13.3cm tall 12.6" wide x 6.7" deep x 5.2" tall



PolyPro units come equipped with:

- Low water detection system
- Saved and stored settings in case of power failure
- Digital three digit display in 0.1°C increments
- Internal temperature uniformity of $\pm 0.4^{\circ}\text{C}$
- Ambient to 99.9°C
- Contain impact resistant see-through polymer covers
- 2 year repair or replacement warranty!