

# Electrophoresis RS - PS - 75

## Power Supply



Revolutionary Science  
Manufacturer of Precision Laboratory Equipment

The Revolutionary Science power supply is the ideal power supply for running routine electrophoresis gels. The RPS-75 is capable of running two gels simultaneously at the same or different voltages. It is equipped with a red warning LED light that illuminates when output current exceeds 350mA. Rapid discharge feature eliminates the possibility of electric shock when the unit is turned off. Two year warranty included.

### RS -PS -75 specifications

Max Power:	52.5 watts
Max Current:	350mA
Output Voltages:	0 volts, 15 volts, 75 volts, 150 volts
Output Jack Sets:	2
Fuse:	Glass encapsulated
Net Weight:	0.68kg (1.5lbs)
Simplicity:	Easy to use slide switches
Dimensions:	7cm tall x 12.7cm wide x 10.2cm deep 2.8" tall x 5" wide x 4" deep



888.322.5722  
www.cplabsafety.com

E-mail: [info@cplabsafety.com](mailto:info@cplabsafety.com)  
Website: [www.calpaclab.com](http://www.calpaclab.com)

ISO 9001  
REGISTERED FIRM