

Pro Water Baths

RS - SP- 100 & RS - PB- 50



Revolutionary Science

Manufacturer of Precision Laboratory Equipment

Quick, durable, compact! The modest sized Mini Pro Bath is perfect for the individual or applications that don't require a large bath. The high impact polypropylene vented cover allows for tall beakers, bottles or flasks to saturate in a precision temperature controlled environment. Two year warranty included.

Mini ProBath RS -PB -50 specifications

Precision:	Internal temperature accuracy of $\pm 2.0^{\circ}\text{C}$
Temperature Range:	Ambient to 99.9°C
Capacity:	0.5L
Net Weight:	1kg (2.3lbs)
Power Requirements:	120 VAC, 50 -60 Hz, 200 W, grounded cord
Durability:	Rugged exterior is resistant to dents and cracks
Portability:	Great for transportation
Dimensions:	
Exterior:	20.3cm diameter x 11.4cm tall 8" diameter x 4.5" tall
Interior:	10.5cm diameter x 6.5cm tall 4.1" diameter x 2.6" tall



With a powerful, precision controlled 1250 watt heating element, the Speed Pro Bath is one of the fastest water baths on the market today. Built from medical grade, hand polished stainless steel, it not only looks sharp, but can also be used for countless applications in the laboratory setting. Includes an unbreakable polycarbonate cover. Two year warranty included.

Speed Pro RS -SP -100 specifications

Precision:	Internal temperature uniformity of $\pm 0.4^{\circ}\text{C}$
Display:	3 digit, 0.1°C increments
Speed:	Temperature reaches boiling in about 10 minutes
Temperature Range:	Ambient to 99.9°C
Capacity:	7.9L
Net Weight:	4kg (9lbs)
Power Requirements:	100-240 VAC, 50 -60 Hz, 1250 W, grounded cord
Durability:	Stainless steel chamber, rugged ABS exterior
Dimensions:	
Exterior:	38cm wide x 30.5cm deep x 25.4cm tall 15" wide x 12" deep x 10" tall
Interior:	23cm diameter x 19cm tall



888.322.5722
www.cplabsafety.com

E-mail: info@cplabsafety.com
Website: www.calpaclab.com

ISO 9001
REGISTERED FIRM