



# TrippNT<sup>TM</sup> MRI

Carts, Bins & Storage



2011 v.1

# TrippNT<sup>™</sup> MRI Catalog

Welcome to TrippNT's New MRI Catalog. It's full of new products designed specifically to make it easier to work, move, store, organize and find supplies. Since many people who work in a MRI facility are on their feet a lot, we've done our homework to design carts, dispensers, bins and holders to eliminate unnecessary steps and effort in your day. Please take time to look through this catalog carefully and ask yourself if it's possible to work easier with some help from one of these revolutionary and affordable products.



## MR Safe

All TrippNT products marked as **MR Safe** are nonconducting, nonmetallic, and nonmagnetic and pose no known hazards in all MR imaging environments.



## MR Conditional

All TrippNT products marked as **MR Conditional** may contain some nonferrous metallic components. These items have been demonstrated to show no known hazards in specified MR environments.

Read more about how we test our MRI compatible products at [TrippNT.com/MRI](https://www.trippnt.com/MRI)



## TrippNT's Triple Trust Guarantee

TrippNT's Triple Trust Guarantee ensures that every MRI Product you choose meets high standards of innovation, testing and standards.



### Innovative Design

Every great product begins with a great idea. Each TrippNT MRI product relies on innovation and original thinking to solve storage and organizing problems elegantly and effectively at the lowest possible price. We promise to creatively apply knowledge, the latest technology and finely-crafted-skill to offer you products that make your day easier.



### Meticulous Testing

Before any new design can wear the TrippNT Triple Trust Shield our Product Testing Team will put it through its paces. Here the product is subjected to rigorous formalized and documented testing procedures to ensure that it will work exactly as intended and will stand up to the abuse that it will face in the field.



### Ridiculously High Standards

Once testing is complete and a design is finalized the new product enters production. From here it will be subject to a documented visual and manual inspection to insure every product meets the high standards of quality specified by TrippNT.

# MRI Core DX™ Lab Cart



The demand for MRI Labs and facilities is growing. Patients must be processed with care and understanding, but also speed. The MRI Core DX is designed to save time by keeping supplies close at hand. Store hearing protection, safety glasses, disposable gloves, protective clothing and more in the three large drawers, six side shelves and bulk storage cabinet of the MRI Core DX. Designed for those who want function and a professional image. Push and pull with comfortable built-in handles.

Product Name	Number	Price
MRI Core DX Lab Cart	51728	\$937.00

# MRI Core DX™ Extended Top




The optional Extended Top MRI Core DX is an extra six inches wide to give more work space and foot room for pushing or pulling the cart.

Product Name	Number	Price
MRI Core DX Extended	51817	\$1,017.00

# MRI Core DX Railed Top

The optional MRI Core DX Railed Top has a 1 3/4" tall rail on all four sides to keep items secure as you move.



Product Name	Number	Price
MRI Core DX Railed Top	50970	\$1042.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call **1.800.874.7768**  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection.  
(Read more on page 3)



# MRI Cayman Lab Island™



The MRI Cayman Lab Island provides an expansive work and storage space with four full feet of movable cabinet. High density polyethylene construction means it will stand up to all types of abuse. With five inch nonferrous rubber wheels, this is one easy-to-move, hard working bench. You'll love the disappearing roll-top door and locking mechanism. It's never in the way.



Product Name	Number	Price
MRI Cayman Lab Island	51821	\$2,230.00
MRI Pad Lock	50971	\$19.00

# MRI Maui Lab Island™



The MRI Maui offers the best of both worlds with liberal storage and work space and a small foot print. Extra large five-inch rubber casters make pushing it a breeze and the high-density polyethylene construction means years of like-new service. 2' x 37 1/2" x 2' WHD



Product Name	Number	Price
MRI Maui Lab Island	51820	\$1,349.00
MRI Pad Lock	50971	\$19.00



## Is your MRI facility ACR certified?

The ACR (American College of Radiology) has a certification program that tests the quality of the images generated by your MRI scanners. Many insurance companies now require their Radiology centers to be ACR certified. Certified centers are allowed to display the ACR seal. Find out more online at: [ACR.org](http://ACR.org)

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection. (Read more on page 3)



# MRI Core SP™ Lab Cart



The MRI Core SP Saves Space by offering a workspace, three large drawers and bulk storage in a 19 1/4 x 15" foot print. The sliding door stays put in any position and saves space since it doesn't swing out. A large bulk storage area holds large items and the clear door lets you see contents from a distance. Smooth rolling casters are great for those who like to store, use and move supplies quickly.

Product Name	Number	Price
MRI Core SP Lab Cart	51726	\$904.00

# MRI Locking Cart



The MRI Locking Cart offers the highest security in the smallest area. With six included pull-out drawers and a unique locking mechanism\* this cart will keep your most valued supplies secure and safe. TrippNT's roll-top door (it disappears when you raise it up) opens smoothly and saves space. Whisk it up and out of the way when you are using supplies and drop it down and lock it when you are away. This cart is great for people who need to protect and secure valuable reagents, tools and equipment.

Product Name	Number	Price
MRI Locking Cart	50729	\$833.00
MRI Pad Lock	50971	\$19.00

\*lock not included

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection.  
(Read more on page 3)





# MRI Extra-Wide Lab Cart



The MRI Extra Wide Lab Cart provides extra-wide mobile work and storage space. Whether you are in a hospital, urgent care facility, physician's office or imaging center safety, professionalism and organization are key to your success.

The MRI Extra Wide Lab Cart provides three oversized drawers and a protected bulk storage area with sliding see-through door and sits atop five inch casters. Its easy mobility and professional appearance will create instant neatness in your MRI Lab or facility.



Product Name	Number	Price
MRI Extra Wide Cart	51823	\$1189.00

# MRI Light Duty Cart



If you need a light-weight cart with heavy-weight function the MRI Light Duty Cart is for you. This super easy to move, ultra-light cart, makes easy work of moving your supplies to your patient's side. This light duty cart has four removable drawers, four side shelves and a bulk storage area.



Product Name	Number	Price
MRI Light Duty Cart	50790	\$515.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call **1.800.874.7768**  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection.  
(Read more on page 3)



# MRI Glove Box Holders

## Top Loading




Product Name	Number	Price
Single	51872	\$36.00
Double	51873	\$46.00
Triple	51874	\$56.00
Triple (side-loading)	51875	\$56.00

*Different Colors available at TrippNT.com*

Mount and identify exam gloves on the wall or bench top with these super-clean yet visual holders for boxes of gloves. The combination of the white PVC body and clear acrylic front gives a hospital-clean appearance with the convenience of a see-through front.

## Side-Loading



Product Name	Number	Price
Single	51876	\$47.00
Double	51877	\$72.00
Triple	51878	\$76.00
Quadruple	51879	\$88.00

Easily exchange glove boxes with these individual compartments that fit nearly any size box of gloves. A large opening makes grabbing gloves convenient. Side-loading means you don't have to remove other boxes of gloves to add or remove new ones.

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



### TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection. (Read more on page 3)




# MRI Storage Boxes



With an open top and wall mountable holes these boxes are truly all purpose on-the-go storage units. They are also roomy and can be piled high with your supplies. The tall back gives extra strength when mounting on the wall although they can sit on the counter, in a refrigerator or cabinet for transport as needed. Install instantly with provided double-sided tape.

Product Name	Number	Price
Small Storage Box	51880	\$39.00
Large Storage Box	51881	\$49.00
<i>Different Colors available at <a href="http://TrippNT.com">TrippNT.com</a></i>		

# MRI Hooks



Plastic hooks hold up to ten pounds on each hook. These MRI Safe hooks install with double-sided foam tape where you want them so your stuff is where you want it. Store patient clothing, purses, scarves, umbrellas, coats, backpacks, towels, cleaning supplies, brooms, anything on a rope, tubing or hangers.

Product Name	Number	Price
Single	51882	\$15.00
Double	51883	\$24.00
Triple	51884	\$32.00
Quadruple	51885	\$39.00
Quintuple	51886	\$45.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection.  
*(Read more on page 3)*





# MRI Storage Bins



51887

Product Name	Number	Price
Eight Fixed Compartment Bin	51887	\$119.00

Sort, store, organize and designate a place for everything with MRI Safe Storage bins. Bins can be mounted on the wall with appropriate hardware or sit freely on a counter top. These bins are made from solid MRI Safe plastic-



51889



51888

10 3/4"

24"

16"

12"

10 1/4"

PVC, are easy to clean and lightweight. Completely durable, strong and tough for years of flexible storage. MRI Bins offer quick inventory assessment.

Product Name	Number	Price
One Shelf and Five Fixed Compartment Bin	51888	\$79.00
One Shelf and Ten Fixed Compartment Bin	51889	\$139.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call **1.800.874.7768**  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection. (Read more on page 3)



# MRI Adjustable Storage Bins



51890

If your supplies ebb and flow you'll find these adjustable bins very versatile. Create your own bin sizes with these adjustable storage bins. Removable dividers make it a snap to adjust bins to fit the amount of supplies on hand.

These bins are made from solid MRI Safe plastic PVC, are easy to clean and lightweight. Completely durable, strong and tough for years of flexible storage.

Product Name	Number	Price
Sixteen Adjustable Compartments	51890	\$109.00



51892

Product Name	Number	Price
Eight Adjustable Compartments and Two Shelves	51891	\$95.00
Fourteen Adjustable Compartments and Two Shelves	51892	\$179.00



51891

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call **1.800.874.7768**  
Free shipping for orders over \$75.




## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection.  
(Read more on page 3)



# MRI Two-Story Bottle Holders



Continuing with TrippNT's innovative and useful products are these new two-story bottle holders. Their 100% plastic construction is MRI Safe, resistant to chemicals, strong and mounts on wall or sits on bench. Installation is instant and they hold 500 mL squeeze bottles, reagents, or chemicals safe and secure. Easy to load and unload without squeezing bottle. Besides the great look these holders sport, they do what they are meant to do well. Hole diameter is 3 1/8" (79 mm).

Product Name	Number	Price
3-Hole Two-Story MRI Bottle Holder	51893	\$28.00
5-Hole Two-Story MRI Bottle Holder	51894	\$38.00
7-Hole Two-Story MRI Bottle Holder	51895	\$48.00
9-Hole Two-Story MRI Bottle Holder	51896	\$58.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection. (Read more on page 3)



# MRI Kleenex Box Holders




51897

This one-box tissue holder attaches instantly with double-faced tape. Just peel off the tape backing and press it in place. Loading the holder only takes seconds by sliding the box in place. The durable, MRI Safe PVC is easy to clean and provides secure access for one-handed tissue dispensing.

Product Name	Number	Price
Industrial Size Holder (100 Count Box)	51897	\$16.00
Regular Size Holder (184 Count Box)	51898	\$17.00
Family Size Holder (280 Count Box)	51899	\$21.00

*More Sizes available at [TrippNT.com](http://TrippNT.com)*

# MRI Paper Towel Holders



Putting paper towels out for use has never been easier or safer. Holder provides easy access while keeping towels dry. Wall mount these holders with included key holes. Overall, this is a secure and simple way to store paper, 'C' or 'Z' fold towels and it's 100% MRI Safe plastic construction.

Product Name	Number	Price
Small Tri-Fold Paper Towel Holder	51900	\$18.00
Medium Tri-Fold Paper Towel Holder	51901	\$28.00
Large Tri-Fold Paper Towel Holder	51902	\$38.00

Each TrippNT MRI product is extensively tested to ensure MRI compliance. (See page 2 for more details)



Order online at [TrippNT.com](http://TrippNT.com) or call 1.800.874.7768  
Free shipping for orders over \$75.



## TrippNT's Triple Trust Guarantee

Each TrippNT product goes through an extensive three step process of design, testing and inspection. (Read more on page 3)





8830 NE 108th Street  
Kansas City, MO 64157



Get a digital copy of this  
catalog online at  
**TrippNT.com/MRI**  
or scan this QR Code  
with your smart phone  
to be taken away into  
the internet.